



j9f6cr3rf

Size selection (Purification) (1X ratio) V.(j9f6cr3rf)



Version 1 is forked from [Size selection \(Purification\)](#).

 Manage mention

Tsu-Chun Hung¹, Yin-Tse Huang¹

¹Kaohsiung Medical University



Yin-Tse Huang

Kaohsiung Medical University



Protocol Info:

Tsu-Chun Hung, Yin-Tse Huang . Size selection (Purification) (1X ratio). **protocols.io**

<https://protocols.io/view/size-selection-purification-1x-ratio-j9f6cr3rf>

Version created by

Yin-Tse Huang

Created: September 13, 2026

Last Modified: September 13, 2026

Protocol Integer ID: 324830

Keywords: purification, size selection, selection

Abstract

Size selection (Purification)



- 1 Calculate the amount of magnetic beads you need and **DISTRIBUTE in a new 1.5 mL eppendorf tube** (Mix the bottle **VIGOROUSLY** to make sure it's homogeneous). For example, if you have 5 samples to clean up, distribute 49.5 μl in a new 1.5 mL eppendorf tube.

	A	B
	Sample	Magnetic beads need
	1	25 μl
	2	50 μl
	3	75 μl
	4	100 μl
	5	125 μl
	6	150 μl
	7	175 μl
	8	200 μl

Note

The ratio is **25 μl of sample** use **25 μl of magnetic beads** to clean up (1 X to keep moderate to long DNA fragments).

Note

Magnetic bead: BeaverBeads™ DNA Select Isolation

- 2 In each  25 μL sample, add  25 μL magnetic beads in.

Note

Pipetting thoroughly every time before adding to make sure all magnetic beads mixed well.

- 3 Mix sample and beads by **gently flicking** the tube for 5mins  00:05:00

5m

Note

Be gentle, mix too aggressive will make DNA fragmented



- 4 Transfer the tube to the 8-strip magnetic rack.
After all the magnetic beads are arrested, remove the supernatant.
- 5 Add  200 μL 80% ethanol, flip the magnetic rack around to clean the pellet.
Wait until all the magnetic beads are arrested, remove the supernatant.

Note

FRESHLY make 80% ethanol

- 6 Repeat the **step 5**.
- 7 Flash spin the tube, and put on the magnetic rack.
Remove **ANY trace of liquid residuals** using 10 μL pipette.

Note

Caution: DO NOT let the beads crack

- 8 Add  15 μL elution buffer.
Mix gently by flicking the tube and make sure all magnetic beads are dissolved in buffer, then flash spin down the sample.
- 9 Incubate the tube in PCR thermal cycler using "**37_incubation**" program for 10min
 00:10:00 .
- 10 Transfer the tube to the magnetic rack.
After all the magnetic beads are arrested, collect the  15 μL supernatant to a 200 μL PCR tube or a 8-strip PCR tube.
- 11 The clean-up sample is now ready for 2' PCR.

10m